

The Mobility Solutions Provider

Scooterland Mobility offers sales, service, repair, hire & spare parts for mobility equipment, products & daily living aids.

Our Workshop Service Comprises:-

Like most elctro-mechanical equipment, mobility products do require servicing to maintain their optimum performance; our recommendation is for servicing to be carried out regularly at yearly intervals for electric mobility scooters, wheelchairs (electric and manual) and powerchairs. At Scooterland we have a fully equipped workshop, with test equipment, spare parts such as tyres, batteries, controllers etc., to undertake servicing and repair work.

Our Service Work Includes?

- ⇒ Carry out a general visual inspection of components
- ⇒ Remove a brush and flush the motor housing assembly of carbon deposit
- ⇒ Check brushes for wear
- ⇒ Remove the wheels and grease the axles(to reduce the risk of a wheel seizing to the axle)
- ⇒ Check wiring connectors to look for corrosion (intermittent faults due to lack of contact)
- ⇒ Check wiring, wiring insulation, for wear or damage (potential malfunction)
- ⇒ Check the charger cord, for wear or damage
- ⇒ Check and adjust tyre pressures
- ⇒ Check front steering for tightness



Our Repair Work includes?

- ⇒ Battery discharge tests
- ⇒ Charger operational tests
- ⇒ Battery replacement
- ⇒ Charger replacement
- ⇒ Tyre or tube replacement
 ⇒ Parts replacement or frame repairs
- ⇒ Motor brushes or brake assembly replacement
- ⇒ Electronic repairs such as throttle potentiometers
- ⇒ Electronic repairs such as speed potentiometers
- ⇒ Troubleshooting fault conditions
- ⇒ Troubleshooting wiring or controller issues
- ⇒ Program diagnostic check
- ⇒ Replacement Controller or Joystick
- ⇒ Controller or Joystick programming




Our After Sales Service

Spare Parts For:-

- ⇒ Wheelchairs (Electric and Manual)
- ⇒ Power and Rehab Chairs
- ⇒ Mobility Scooters
- ⇒ Mobility Aids (Walkers)
- ⇒ Batteries and Chargers
- ⇒ Electronics Controllers & Joysticks

Accessories:-

- ⇒ Rear Bags
- ⇒ Canopies
- ⇒ Walking Frames





SALES * SERVICE * REPAIR * HIRE

07 5539 2141



sales@scooterland.com.au

www.scooterland.com.au